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
- Albert MS, et al. **The diagnosis of mild cognitive impairment due to Alzheimer's disease: recommendations from the National Institute on Aging-Alzheimer's Association workgroups on diagnostic guidelines for Alzheimer's disease.** *Alzheimers Dement* 7(3):270-9 (2011). <https://doi.org/10.1016/j.jalz.2011.03.008>
- DGN + DGLN. **S1-Leitlinie „Lumbalpunktion und Liquordiagnostik“.** AWMF-Register-Nr. 030/141 (2019). <https://www.awmf.org/leitlinien/detail/II/030-141.html>
- DGPPN + DGN. **S3-Leitlinie „Demenzen“.** AWMF-Register-Nr. 038-013 (2016). <https://www.awmf.org/leitlinien/detail/II/038-013.html>
- Dubois B, et al. **Clinical diagnosis of Alzheimer's disease: recommendations of the International Working Group.** *Lancet Neurol* S1474-4422(21)00066-1 (2021). [https://doi.org/10.1016/S1474-4422\(21\)00066-1](https://doi.org/10.1016/S1474-4422(21)00066-1)
- Jack Jr CR, et al. **Introduction to revised criteria for the diagnosis of Alzheimer's disease: National Institute on Aging and the Alzheimer Association Workgroups.** *Alzheimers Dement* 7(3):257-62 (2011). <https://doi.org/10.1016/j.jalz.2011.03.004>
- McKhann GM, et al. **The diagnosis of dementia due to Alzheimer's disease: Recommendations from the National Institute on Aging-Alzheimer's Association workgroups on diagnostic guidelines for Alzheimer's disease.** *Alzheimers Dement* 7(3):263-9 (2011). <https://doi.org/10.1016/j.jalz.2011.03.005>
- Sorbi S, et al. **EFNS-ENS Guidelines on the diagnosis and management of disorders associated with dementia.** *Eur J Neurol* 19(9):1159-79 (2012). <https://doi.org/10.1111/j.1468-1331.2012.03784.x>
- Sperling RA, et al. **Toward defining the preclinical stages of Alzheimer's disease: Recommendations from the National Institute on Aging-Alzheimer's Association workgroups on diagnostic guidelines for Alzheimer's disease.** *Alzheimers Dement* 7(3):280-92 (2011). <http://doi.org/10.1016/j.jalz.2011.03.003>

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Lumbalpunktion

- Babapour Mofrad R, et al. **Cerebrospinal fluid collection: an informative animation video for patients and caregivers.** *Alzheimers Dement (Amst)* 11(1):435-438 (2019). <https://doi.org/10.1016/j.dadm.2019.04.005>
- Babapour Mofrad R, et al. **Lumbar puncture in patients with neurologic conditions.** *Alzheimers Dement (Amst)* 8(1):108-110 (2017). <https://doi.org/10.1016/j.dadm.2017.04.008>
- DGN + DGLN. **S1-Leitlinie „Lumbalpunktion und Liquordiagnostik“.** AWMF-Register-Nr. 030/141 (2019). <https://www.awmf.org/leitlinien/detail/II/030-141.html>
- Engelborghs S, et al. **Consensus guidelines for lumbar puncture in patients with neurological diseases.** *Alzheimers Dement (Amst)* 8(1):111-126 (2017). <https://doi.org/10.1016/j.dadm.2017.04.007>
- Rembach A, et al. **Alzheimer's disease cerebrospinal fluid biomarkers are not influenced by gravity drip or aspiration extraction methodology.** *Alzheimers Res Ther* 7(1):71 (2015). <https://doi.org/10.1186/s13195-015-0157-7>
- Shaw LM, et al. **Appropriate use criteria for lumbar puncture and cerebrospinal fluid testing in the diagnosis of Alzheimer's disease.** *Alzheimers Dement* 14(11):1505-1521 (2018). <https://doi.org/10.1016/j.jalz.2018.07.220>
- Teunissen CE, et al. **A consensus protocol for the standardization of cerebrospinal fluid collection and biobanking.** *Neurology.* 73(22):1914-22 (2009). <https://doi.org/10.1212/WNL.0b013e3181c47cc2>

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Probenhandhabung (SOP)

- del Campo M, et al. **Recommendations to standardize preanalytical confounding factors in Alzheimer's and Parkinson's disease cerebrospinal fluid biomarkers: an update** *Biomark Biomark Med* 6(4):419-30 (2012). <https://doi.org/10.2217/bmm.12.46>
- Hansson O, et al. **The impact of preanalytical variables on measuring cerebrospinal fluid biomarkers for Alzheimer's disease diagnosis: a review.** *Alzheimers Dement* 14(10):1313-1333 (2018). <https://doi.org/10.1016/j.jalz.2018.05.008>
- Hansson O, et al. **Pre-analytical protocol for measuring Alzheimer's disease biomarkers in fresh CSF.** *Alzheimers Dement (Amst)* 12(1):e12137 (2020). <https://doi.org/10.1002/dad2.12137> 
- Hansson O, et al. **The Alzheimer's Association international guidelines for handling of cerebrospinal fluid for routine clinical measurements of amyloid β and tau.** *Alzheimers Dement* 17(9):1575-1582 (2021). <https://doi.org/10.1002/alz.12316> 
- Vanderstichele H, et al. **Standardization of preanalytical aspects of cerebrospinal fluid biomarker testing for Alzheimer's disease diagnosis: A consensus paper from the Alzheimer's Biomarkers Standardization Initiative.** *Alzheimers Dement* 8(1):65-73 (2012). <https://doi.org/10.1016/j.jalz.2011.07.004>

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Präanalytik

- Andreasen N, et al. **Cerebrospinal Fluid β -Amyloid(1-42) in Alzheimer Disease. Differences Between Early- and Late-Onset Alzheimer Disease and Stability During the Course of Disease.** Arch Neurol 56(6):673-80 (1999). <https://doi.org/10.1001/archneur.56.6.673>
- Bellomo G, et al. **Measurement of CSF core Alzheimer disease biomarkers for routine clinical diagnosis: do fresh vs frozen samples differ?** Alzheimers Res Ther 12:121 (2020). <https://doi.org/10.1186/s13195-020-00689-0>
- Bjerke M, et al. **Confounding Factors Influencing Amyloid Beta Concentration in Cerebrospinal Fluid.** Int J Alzheimers Dis 2010:986310 (2010). <https://doi.org/10.4061/2010/986310>
- Cicognola C, et al. **No diurnal variation of classical and candidate biomarkers of Alzheimer's disease in CSF.** Mol Neurodegener 11(1):65 (2016). <https://doi.org/10.1186/s13024-016-0130-3>
- del Campo M, et al. **Recommendations to standardize preanalytical confounding factors in Alzheimer's and Parkinson's disease cerebrospinal fluid biomarkers: an update.** Biomark Med 6(4):419-30 (2012). <https://doi.org/10.2217/bmm.12.46>
- Gervaise-Henry C, et al. **Cerebrospinal Fluid $A\beta_{42}/A\beta_{40}$ as a Means to Limiting Tube- and Storage-Dependent Pre-Analytical Variability in Clinical Setting.** J Alzheimers Dis 57(2):437-445 (2017). <https://doi.org/10.3233/JAD-160865>
- Hansson O, et al. **The impact of preanalytical variables on measuring cerebrospinal fluid biomarkers for Alzheimer's disease diagnosis: a review.** Alzheimers Dement 14(10):1313-1333 (2018). <https://doi.org/10.1016/j.jalz.2018.05.008>
- Hansson O, et al. **Advantages and disadvantages of the use of the CSF Amyloid β ($A\beta$) 42/40 ratio in the diagnosis of Alzheimer's Disease.** Alzheimers Res Ther 11(1):34 (2019). <https://doi.org/10.1186/s13195-019-0485-0>

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Präanalytik

- Hansson O, et al. **Pre-analytical protocol for measuring Alzheimer's disease biomarkers in fresh CSF.** *Alzheimers Dement (Amst)* 12(1):e12137 (2020). <https://doi.org/10.1002/dad2.12137>
- Hansson O, et al. **The Alzheimer's Association international guidelines for handling of cerebrospinal fluid for routine clinical measurements of amyloid β and tau.** *Alzheimers Dement* 17(9):1575-1582 (2021). <https://doi.org/10.1002/alz.12316>
- Janelidze S, et al. **CSF A β 42/A β 40 and A β 42/A β 38 ratios: better diagnostic markers of Alzheimer disease.** *Ann Clin Transl Neurol* 3(3):154-65 (2016). <https://doi.org/10.1002/acn3.274>
- Janelidze S, et al. **Towards a unified protocol for handling of CSF before β -amyloid measurements.** *Alzheimers Res Ther* 11(1):63 (2019). <https://doi.org/10.1186/s13195-019-0517-9>
- Kuhlmann J, et al. **CSF A β ₁₋₄₂ – an excellent but complicated Alzheimer's biomarker – a route to standardisation.** *Clin Chim Acta* 467:27-33 (2017). <https://doi.org/10.1016/j.cca.2016.05.014>
- Le Bastard N, et al. **Importance and Impact of Preanalytical Variables on Alzheimer Disease Biomarker Concentrations in Cerebrospinal Fluid.** *Clin Chem* 61(5):734-43 (2015). <https://doi.org/10.1373/clinchem.2014.236679>
- Lehmann S, et al. **Impact of harmonization of collection tubes on Alzheimer's disease diagnosis.** *Alzheimers Dement* 10(5 Suppl):S390-S394.e2 (2014). <https://doi.org/10.1016/j.jalz.2013.06.008>
- Leitão MJ, et al. **Chasing the effects of pre-analytical confounders – a multicenter study on CSF-AD biomarkers.** *Front Neurol* 6:153 (2015). <https://doi.org/10.3389/fneur.2015.00153>

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Präanalytik

- Lewczuk P, et al. **Effect of Sample Collection Tubes on Cerebrospinal Fluid Concentrations of Tau Proteins and Amyloid β Peptides.** Clin Chem 52(2):332-4 (2006). <https://doi.org/10.1373/clinchem.2005.058776>
- Perret-Liaudet A, et al. **Risk of Alzheimer's Disease Biological Misdiagnosis Linked to Cerebrospinal Collection Tubes.** J Alzheimers Dis 31(1):13-20 (2012). <https://doi.org/10.3233/JAD-2012-120361>
- Rembach A, et al. **Alzheimer's disease cerebrospinal fluid biomarkers are not influenced by gravity drip or aspiration extraction methodology.** Alzheimers Res Ther 7(1):71 (2015). <https://doi.org/10.1186/s13195-015-0157-7>
- Schauer SP, et al. **Preanalytical approaches to improve recovery of amyloid- β peptides from CSF as measured by immunological or mass spectrometry-based assays.** Alzheimers Res Ther 10(1):118 (2018). <https://doi.org/10.1186/s13195-018-0445-0>
- Schoonenboom NS, et al. **Effects of Processing and Storage Conditions on Amyloid β (1–42) and Tau Concentrations in Cerebrospinal Fluid: Implications for Use in Clinical Practice.** Clin Chem 51(1):189-95 (2005). <https://doi.org/10.1373/clinchem.2004.039735>
- Toombs J, et al. **Identification of an important potential confound in CSF AD studies: aliquot volume.** Clin Chem Lab Med 51(12):2311-7 (2013). <https://doi.org/10.1515/cclm-2013-0293>
- Toombs J, et al. **Amyloid β peptides are differentially vulnerable to preanalytical surface exposure, an effect incompletely mitigated by the use of ratios.** Alzheimers Dement (Amst) 10:311–321 (2018). <https://doi.org/10.1016/j.dadm.2018.02.005>
- Vanderstichele H, et al. **Standardization of preanalytical aspects of cerebrospinal fluid biomarker testing for Alzheimer's disease diagnosis: A consensus paper from the Alzheimer's Biomarkers Standardization Initiative.** Alzheimers Dement 8(1):65-73 (2012). <https://doi.org/10.1016/j.jalz.2011.07.004>

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Präanalytik

- Vanderstichele H, et al. **Optimized Standard Operating Procedures for the Analysis of Cerebrospinal Fluid A β 42 and the Ratios of A β Isoforms Using Low Protein Binding Tubes.** J Alzheimers Dis 53(3):1121-32 (2016). <https://doi.org/10.3233/JAD-160286>
- Vanderstichele H, et al. **Recommendations for cerebrospinal fluid collection for the analysis by ELISA of neurogranin trunc P75, α -synuclein, and total tau in combination with A β (1-42)/A β (1-40).** Alzheimers Res Ther 9(1):40 (2017). <https://doi.org/10.1186/s13195-017-0265-7>
- Vanderstichele H, et al. **Standardization of measurement of β -amyloid₍₁₋₄₂₎ in cerebrospinal fluid and plasma.** Amyloid 7(4):245-58 (2000). <https://doi.org/10.3109/13506120009146438>
- Willemse E, et al. **How to handle adsorption of cerebrospinal fluid amyloid β (1-42) in laboratory practice? Identifying problematic handlings and resolving the issue by use of the A β 42/A β 40 ratio.** Alzheimers Dement 13(8):885-892 (2017). <http://doi.org/10.1016/j.jalz.2017.01.010>
- Zimmermann R, et al. **Preanalytical Sample Handling and Sample Stability Testing for the Neurochemical Dementia Diagnostics.** J Alzheimers Dis 25(4):739-45 (2011). <https://doi.org/10.3233/JAD-2011-110212>

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- Joint Research Centre. <https://crm.jrc.ec.europa.eu/?q=amyloid>
- Andreasson U, et al. **Commutability of the certified reference materials for the standardization of β -amyloid 1-42 assay in human cerebrospinal fluid: lessons for tau and β -amyloid 1-40 measurements.** Clin Chem Lab Med 56(12):2058-2066 (2018). <https://doi.org/10.1515/cclm-2018-0147>
- Bjerke M, et al. **Assessing the commutability of reference material formats for the harmonization of amyloid- β measurements.** Clin Chem Lab Med 54(7):1177-91 (2016). <https://doi.org/10.1515/cclm-2015-0733>
- Boulo S, et al. **First amyloid β 1-42 certified reference material for re-calibrating commercial immunoassays.** Alzheimers Dement 16(11):1493-1503 (2020). <https://doi.org/10.1002/alz.12145>
- Kuhlmann J, et al. **CSF $A\beta_{1-42}$ – an excellent but complicated Alzheimer's biomarker – a route to standardisation.** Clin Chim Acta 467:27-33 (2017). <https://doi.org/10.1016/j.cca.2016.05.014>

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Interpretationsalgorithmen

- Dubois B, et al. **Clinical diagnosis of Alzheimer's disease: recommendations of the International Working Group.** Lancet Neurol S1474-4422(21)00066-1 (2021). [https://doi.org/10.1016/S1474-4422\(21\)00066-1](https://doi.org/10.1016/S1474-4422(21)00066-1)
- Ekman U, et al. **The A/T/N biomarker scheme and patterns of brain atrophy assessed in mild cognitive impairment.** Sci Rep 8(1):8431 (2018). <https://doi.org/10.1038/s41598-018-26151-8>
- Jack Jr CR, et al. **A/T/N: An unbiased descriptive classification scheme for Alzheimer disease biomarkers.** Neurology 87(5):539-47 (2016). <https://doi.org/10.1212/WNL.0000000000002923>
- Jack Jr CR, et al. **NIA-AA Research Framework: Toward a biological definition of Alzheimer's disease.** Alzheimers Dement 14(4):535-562 (2018). <https://doi.org/10.1016/j.jalz.2018.02.018>
- Lehmann S, et al. **A diagnostic scale for Alzheimer's disease based on cerebrospinal fluid biomarker profiles.** Alzheimers Res Ther 6(3):38 (2014). <https://doi.org/10.1186/alzrt267>
- Lehmann S, et al. **Relevance of A β 42/40 Ratio for Detection of Alzheimer Disease Pathology in Clinical Routine: The PLM_R Scale.** Front Aging Neurosci 10:138 (2018). <https://doi.org/10.3389/fnagi.2018.00138>
- Lewczuk P, et al. **Neurochemical dementia diagnostics: a simple algorithm for interpretation of the CSF biomarkers.** J Neural Transm (Vienna) 116(9):1163-7 (2009). <https://doi.org/10.1007/s00702-009-0277-y>
- Lewczuk P, et al. **Validation of the Erlangen Score Algorithm for the Prediction of the Development of Dementia due to Alzheimer's Disease in Pre-Dementia Subjects.** J Alzheimers Dis 48(2):433-41 (2015). <https://doi.org/10.3233/JAD-150342>
- Delaby C, et al. **Clinical reporting following the quantification of cerebrospinal fluid biomarkers in Alzheimer's disease: An international overview.** Alzheimers Dement (2021). <https://doi.org/10.1002/alz.12545>

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